

Peristeen® Plus transanal irrigation with cone catheter

New!



Peristeen® Plus

Introducing Peristeen® Plus transanal irrigation (TAI) with cone catheter

At Coloplast, we aim to live up to the expectations you and your patients have. And so, we have introduced a cone catheter to our industry leading transanal irrigation system, increasing the options available to patients with bowel dysfunction.

Read on to find out more about the new Peristeen Plus cone catheter and how it will benefit users.

The Peristeen Plus cone catheter has been developed with input from colorectal surgeons and nurses regarding shape, curvature, length and flexibility

Rounded tip

The rounded tip with two large water eyelets is designed to minimise discomfort during irrigation



Hydrophilic coating

The hydrophilic coating makes the cone easy to insert

Dual holding positions

The dual holding positions enable users to use either the finger grip or cone base to insert the cone



Designed with extra space

Designed with extra space for grip/fingers under the cone base



Peristeen Plus TAI with cone catheter order codes

Order code	Description	Contents
29160	Peristeen Plus system with cone catheter (incl. toiletry bag)	1 control unit, 2 cone catheters, 1 water bag (incl. screw top), 2 straps, 1 tube
29161	Peristeen Plus system with cone catheter (excl. toiletry bag)	1 control unit, 2 cone catheters, 1 water bag (incl. screw top), 2 straps, 1 tube
29162	Peristeen Plus cone catheter accessory unit	15 cone catheters, 1 water bag (excl. screw top)

To find out more about Peristeen Plus, please call 0800 374 654 or email hcp@coloplastcharter.co.uk



CV1074N



Ostomy Care / Continence Care / Wound & Skin Care / Interventional Urology

Coloplast Ltd, Nene Hall, Peterborough Business Park, Peterborough PE2 6FX
www.coloplast.co.uk The Coloplast logo is a registered trademark of Coloplast A/S. © 2022-03. All rights reserved Coloplast A/S. PM-20955

